6003 Chapel Hill Road · Suite # 153 · Raleigh, NC 27607

Certificate of Analysis

Product: Purified Human IgM, Whole Molecule

Part Number: Hu-006-M Lot Number: 68-159-110519

Concentration: $8.80 \text{ mg/ml} (E^{1\%} \text{ at } 280 \text{ nm} = 11.8)$

Amount: 1.0 mg

Form: Clear, colorless liquid Purity: > 95 % by SDS PAGE

Host: Human

Source: Human Myeloma Plasma

Buffer: 50 mM TRIS, pH 8.0, 0.2 M Sodium Chloride

Preservative: 0.05 % (w/v) Sodium Azide

Storage: 2-8 °C

Expiration Date: November 5, 2022

Specificity: Based on IEP, precipitin bands occur when reacted against antibodies to:

• Whole Human serum

• Human IgM

Kappa light chain

Cross Reactivity: Based on IEP, this protein does not form precipitin bands against antibodies

to:

• Human IgA, IgG, IgD, IgE or lambda light chain

Note: To be used as a blocking reagent, normal control, or standard in most

immunoassay protocols. The optimal working dilution should be determined

by the investigator.

Disclaimer: For in vitro Laboratory Use Only. Not for diagnostic or therapeutic use. Not

for human or animal consumption. The proper selection and use of our products is the sole responsibility of the end user, and therefore we can offer no guarantee to a specific experimental outcome. Suggested applications of our products are not recommendations to use our products in violation of any

patent or as a license under any patent of ImmunoReagents, Inc.

Product may not be resold or modified for resale without prior written

approval of ImmunoReagents, Inc.

Phone: 919-831-2240 · www.ImmunoReagents.com · Fax: 919-831-2240

6003 Chapel Hill Road · Suite # 153 · Raleigh, NC 27607

Testing Results

PRODUCT: Purified Human IgM, Whole Molecule

PART NUMBER: Hu-006-M **LOT NUMBER:** 68-159-110519

<u>TESTS</u> <u>SPECIFICATIONS</u> <u>RESULTS</u>

Appearance: Clear, colorless liquid PASS

Purity: Two predominant bands at 70 kDa and 24 kDa

based on SDS-PAGE PASS

SDS-PAGE

1 2

55kDa

31kDa

21kDa

- 1) Hu-006-M, Lot # 68-159-110519
- 2) Mol Weight Standard

Date: November 8, 2019 Signature: Quality Control

Phone: 919-831-2240 • www.ImmunoReagents.com • Fax: 919-831-2240