6003 Chapel Hill Road Suite 153 Raleigh, NC 27607

## **Product Specifications**

Product: Chicken anti-Mouse IgG (H&L) - Affinity Pure

Description: Chicken anti-Mouse IgG (H&L) - Affinity Pure, TRITC Conjugate, min x

w/human or rabbit lgG/serum proteins

Part Number: CkxMu-003-ERHOX

Concentration: 1.0 mg/ml (E 1% at 280 nm = 13.2)

Amount: 0.5 mg

Conjugate: Tetramethylrhodamine-5-isothiocyanate (TRITC) Amax = 550 nm; Emax =

570 nm

Form: Lyophilized

Purification: Affinity purified using solid phase Mouse IgG

Purity: Affinity purified antibody is > 95% based on SDS-PAGE

Host: Chicken

Immunogen: Purified Mouse IgG, whole molecule

Buffer: 10 mM Sodium Phosphate, 0.15 M Sodium Chloride, pH 7.2, 1 % (w/v)

BSA, Protease/IgG free

Preservative: 0.05% (w/v) Sodium Azide

Reconstitution: Rehydrate with 0.55 ml of deionized water and let stand 30 minutes at

room temperature to dissolve. (Product has been overfilled to ensure

complete recovery.) Centrifuge to remove any particulates. Prepare fresh

working dilution daily.

Storage: 2-8 °C

Shelf Life: 1 year from date of receipt. Prepare working dilution prior to use and then

discard.

Specificity: Based on IEP, this antibody reacts with:

heavy (γ) chains on mouse IgG

light chains on all mouse immunoglobulins

Cross Reactivity: Based on IEP, no reactivity is observed to:

non-immunoglobulin mouse serum immunoglobulins

IgG/serum from human or rabbit

Country of Origin: Chicken Ig fraction was prepared using eggs from healthy hens of US

origin.

Applications: Flow Cytometry

Immunofluorescence

Phone: 919-831-2240 Fax: 919-831-2240 email: info@immunoreagents.com

**ELISA** 

Immunomicroscopy

Disclaimer:

For *in vitro* Laboratory Use Only. Not for diagnostic or therapeutic use. Not for human or animal consumption. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license under any patent of ImmunoReagents, Inc. **Product may not be resold or modified for resale without prior written approval of ImmunoReagents, Inc.** 

Phone: 919-831-2240 Fax: 919-831-2240 email: info@immunoreagents.com