6003 Chapel Hill Road Suite 153 Raleigh, NC 27607

Product Specifications

Product: Goat anti-Guinea Pig IgG (H&L) - Affinity Pure

Description: Goat anti-Guinea Pig IgG (H&L) - Affinity Pure, DyLight®488 Conjugate,

min w/bovine, chicken, goat, hamster, horse, human, mouse, rabbit, rat,

sheep Serum

Part Number: GtxGp-003-F488NHSX

Concentration: 1.0 mg/ml (E 1% at 280 nm = 13.0)

Amount: 1.0 mg

Conjugate: DyLight® 488 (Ex = 493 nm; Em = 518 nm)

Form: Lyophilized

Purification: Affinity purified using solid phase Guinea Pig IgG

Purity: Affinity purified antibody is > 95% based on SDS-PAGE

Host: Goat

Immunogen: Purified Guinea Pig IgG, whole molecule

Buffer: 10 mM Sodium Phosphate, 0.15 M Sodium Chloride, pH 7.2, 1 % (w/v)

BSA, Protease/IgG free

Preservative: 0.05% (w/v) Sodium Azide

Reconstitution: Rehydrate with 1.1 ml of deionized water and let stand 30 minutes at room

temperature to dissolve. (Product has been overfilled to ensure complete recovery.) Centrifuge to remove any particulates. Prepare fresh working

dilution daily.

Storage: Store freeze-dried powder at 2-8 °C.

Shelf Life: Product is stable for up to 4 weeks at 2-8 °C after rehydration. For

extended storage after rehydration, add an equal volume of glycerol and

store at -20°C.

Specificity: Based on IEP, this antibody reacts with:

heavy (γ) chains on guinea pig IgG

light chains on all guinea pig immunoglobulins

Cross Reactivity: Based on IEP, no reactivity is observed to:

non-immunoglobulin guinea pig serum immunoglobulins

• IgG from bovine, chicken, goat, hamster, horse, human, mouse,

rabbit, rat or sheep

Country of Origin: Goat serum was obtained from healthy animals of US origin and under the

care of a registered veterinarian.

Phone: 919-831-2240 Fax: 919-831-2240 email: info@immunoreagents.com

Applications: Flow Cytometry

Immunofluorescence

ELISA

Immunomicroscopy

Disclaimer: For in vitro Laboratory Use Only. Not for diagnostic or therapeutic use. Not

for human or animal consumption. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license under any patent of ImmunoReagents, Inc. **Product may not**

be resold or modified for resale without prior written approval of

ImmunoReagents, Inc.

Trademark: DyLight® is a trademark of Thermo Fisher Scientific, Inc. and its

subsidiaries.

Phone: 919-831-2240 Fax: 919-831-2240 email: info@immunoreagents.com