6003 Chapel Hill Road Suite 153 Raleigh, NC 27607

Product Specifications

Product: Goat anti-Mouse IgG (H&L), F(ab)'2 Fragment - Affinity Pure

Description: Goat anti-Mouse IgG (H&L), F(ab)'2 Fragment - Affinity Pure, HRP

Conjugate

Part Number: GtxMu-003-GHRPX

Concentration: 1.0 mg/ml (E 1% at 280 nm = 13.0)

Amount: 0.5 mg

Conjugate: Horseradish Peroxidase

Form: Lyophilized

Purification: Affinity purified using solid phase Mouse IgG

Purity: > 90% based on SDS-PAGE Small amounts of intact IgG may be present.

Host: Goat

Immunogen: Purified Mouse IgG, whole molecule

Buffer: 10 mM Sodium Phosphate, 0.15 M Sodium Chloride, pH 7.2, 1 % (w/v)

BSA, Protease/IgG free

Preservative: 0.1% (v/v) Kathon CG

Reconstitution: Rehydrate with 0.55 ml of deionized water and let stand 30 minutes at

room temperature to dissolve. (Product has been overfilled to ensure complete recovery.) Centrifuge to remove any particulates. Prepare fresh

working dilution daily.

Storage: Store freeze-dried powder at 2-8 °C.

Shelf Life: Store lyophilized material at 2-8 °C. For long term storage after

reconstitution, dilute with 50% glycerol and store at -20 °C as a liquid.

Specificity: Based on IEP, this antibody reacts with:

heavy (γ) chains on mouse IgG

• light chains on all mouse immunoglobulins

Cross Reactivity: Based on IEP, no reactivity is observed to:

non-immunoglobulin mouse serum immunoglobulins

Country of Origin: Goat serum was obtained from healthy animals of US origin and under the

care of a registered veterinarian.

Applications: Western Blot

ELISA

Immunohistochemistry (IHC)

Disclaimer:

For in vitro Laboratory Use Only. Not for diagnostic or therapeutic use. Not

Phone: 919-831-2240 Fax: 919-831-2240 email: info@immunoreagents.com

for human or animal consumption. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license under any patent of ImmunoReagents, Inc. **Product may not be resold or modified for resale without prior written approval of ImmunoReagents, Inc.**

Phone: 919-831-2240 Fax: 919-831-2240 email: info@immunoreagents.com