6003 Chapel Hill Road Suite 153 Raleigh, NC 27607

## **Product Specifications**

Product: Goat anti-Rat lgG (H&L) - Affinity Pure

Description: Goat anti-Rat IgG (H&L) - Affinity Pure, HRP Conjugate, min x w/human or

mouse IgG

Part Number: GtxRt-003-EHRPX

Concentration: 1.0 mg/ml (E 1% at 280 nm = 13.0)

Amount: 1.0 mg

Conjugate: Horseradish Peroxidase

Form: Lyophilized

Purification: Affinity purified using solid phase Rat IgG

Purity: Affinity purified antibody is > 95% based on SDS-PAGE

Host: Goat

Immunogen: Purified Rat IgG, whole molecule

Buffer: 10 mM Sodium Phosphate, 0.15 M Sodium Chloride, pH 7.2, 1 % (w/v)

BSA, Protease/IgG free

Preservative: 0.1% (v/v) Kathon CG

Reconstitution: Rehydrate with 1.1 ml of deionized water and let stand 30 minutes at room

temperature to dissolve. (Product has been overfilled to ensure complete recovery.) Centrifuge to remove any particulates. Prepare fresh working

dilution daily.

Storage: Store freeze-dried powder at 2-8 °C.

Shelf Life: Store lyophilized material at 2-8 °C. For long term storage after

reconstitution, dilute with 50% glycerol and store at -20 °C as a liquid.

Specificity: Based on IEP, this antibody reacts with:

heavy (γ) chains on rat IgG

light chains on all rat immunoglobulins

Cross Reactivity: Based on IEP, no reactivity is observed to:

non-immunoglobulin rat serum proteins

IgG from human or mouse

Country of Origin: Goat serum was obtained from healthy animals of US origin and under the

care of a registered veterinarian.

Applications: Western Blot

**ELISA** 

Immunohistochemistry (IHC)

Phone: 919-831-2240 Fax: 919-831-2240 email: info@immunoreagents.com

Disclaimer:

For *in vitro* Laboratory Use Only. Not for diagnostic or therapeutic use. Not for human or animal consumption. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license under any patent of ImmunoReagents, Inc. **Product may not be resold or modified for resale without prior written approval of ImmunoReagents, Inc.** 

Phone: 919-831-2240 Fax: 919-831-2240 email: info@immunoreagents.com