

6003 Chapel Hill Road · Suite # 153 · Raleigh, NC 27607

Certificate of Analysis

Product:	Rabbit anti-Mouse IgG (H&L), Affinity Pure	
Part Number: Lot Number: Concentration: Amount: Form: Purification: Purity: Host: Immunogen: Buffer: Preservative: Storage: Expiration Date: Specificity:	RbxMu-003-D 61-152-043020 $10.0 \text{ mg/ml} (E^{1\%} \text{ at } 280 \text{ nm} = 13.0)$ 2.0 mg Clear, colorless liquid, $0.2 \mu \text{m}$ filtered Affinity purified using solid phase Mouse IgG $\geq 95\%$ based on SDS-PAGE Rabbit Purified Mouse IgG, whole molecule 10 mM Sodium Phosphate, $0.15 M$ Sodium Chloride, pH 7.2 0.05 % (w/v) Sodium Azide $2-8 \degree \text{C}$ April 30, 2023 Based on IEP, this antibody reacts with: • heavy (γ) chains on mouse IgG • light chains on all mouse immunoglobulins	
Cross Reactivity:	Based on IEP, no reactivity is observed to:non-immunoglobulin mouse serum proteins	
Country of Origin:	Rabbit serum was obtained from healthy animals of US origin and under the care of a registered veterinarian.	
Applications:	This antibody is suitable for use in most immunoassay formats for the detection of Mouse IgG or as a capture antibody in ELISA type formats. This antibody is suitable for use in Western ImmunoBlots and ImmunoHistoChemistry (IHC). This antibody will cross react with IgG from other species. The optimal working dilution should be determined by the investigator.	
Disclaimer:	<i>For in vitro Laboratory Use Only.</i> Not for diagnostic or therapeutic use. Not for mouse or animal consumption. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license under any patent of ImmunoReagents, Inc.	
	Product may not be resold or modified for resale without prior written approval of ImmunoReagents, Inc.	

*ImmunoReagents,Inc.

6003 Chapel Hill Road · Suite # 153 · Raleigh, NC 27607

Testing Results

PRODUCT: PART NUMBER: LOT NUMBER:	Rabbit anti-Mouse IgG (H&L), Affinity Pure RbxMu-003-D 61-152-043020	
<u>TESTS</u>	SPECIFICATIONS	<u>RESULTS</u>
Appearance:	Clear and Colorless	PASS
Purity:	Two predominant bands at 50 and 24 kDa based on SDS PAGE	PASS
Specificity:	 Based on IEP, this antibody has precipitin bands with: heavy (γ) chains on mouse IgG light chains on all mouse immunoglobulins 	PASS PASS
Cross Reactivity:	 Based on IEP, this antibody does not have precipitin bands to: non-immunoglobulin mouse serum proteins 	PASS
SDS PAGE 1 2 55k 31k 21k	Da Normal Mouse Serum	
1) RbxMu-00	3-D, Lot # 61-152-043020	

Mol Weight Standard

Date: May 6, 2020

Signature:

Quality Control