6003 Chapel Hill Road · Suite # 153 · Raleigh, NC 27607

Certificate of Analysis

Product: Goat anti-Human IgA (α chain), Affinity Pure

Part Number: GtxHu-001-D Lot Number: 33-66-112613

Concentration: 9.49 mg/ml (E $^{1\%}$ at 280 nm = 13.0)

Amount: 2.0 mg

Form: Clear, colorless liquid, 0.2 µm filtered

Purification: Affinity purified using solid phase Human IgA

Purity: > 95 % based on SDS-PAGE

Host: Goat

Immunogen: Purified Human IgA, α chain

Buffer: 10 mM Sodium Phosphate, 0.15 M Sodium Chloride, pH 7.2

Preservative: 0.05 % (w/v) Sodium Azide

Storage: 2-8 °C

Expiration Date: June 22, 2024

Specificity: Based on IEP, this antibody reacts with:

• heavy (α) chains on Human IgA

Cross Reactivity: Based on IEP, no reactivity is observed to:

light chains on all Human immunoglobulinsnon-immunoglobulin Human serum proteins

Country of Origin: Goat serum was obtained from healthy animals of US origin and under the

care of a registered veterinarian.

Applications: Suitable for use in: ELISA, WB, IHC, ICC, CL, FACS, IM, IF

The optimal working dilution should be determined by the investigator.

Disclaimer: For in vitro Laboratory Use Only. Not for diagnostic or therapeutic use. Not

for Human or animal consumption. The proper selection and use of our products is the sole responsibility of the end user, and therefore we can offer no guarantee to a specific experimental outcome. Suggested applications of our products are not recommendations to use our products in violation of any

patent or as a license under any patent of ImmunoReagents, Inc.

Product may not be resold or modified for resale without prior written

approval of ImmunoReagents, Inc.

Phone: 919-831-2240 · www.ImmunoReagents.com · Fax: 919-831-2240

6003 Chapel Hill Road · Suite # 153 · Raleigh, NC 27607

Testing Results

PRODUCT: Goat anti-Human IgA (α chain), Affinity Pure

PART NUMBER: GtxHu-001-D **LOT NUMBER:** 33-66-112613

<u>TESTS</u> <u>SPECIFICATIONS</u> <u>RESULTS</u>

Appearance: Clear, colorless liquid PASS

Purity: Two predominant bands at 50 and 24 kDa based on SDS-PAGE PASS

Specificity: Based on IEP, this antibody has precipitin bands with:


heavy (α) chains on Human IgAPASS

Cross Reactivity: Based on IEP, this antibody does not have precipitin bands to:

light chains on all Human immunoglobulins
non-immunoglobulin Human serum proteins
PASS

SDS-PAGE

Immunoelectrophoresis (IEP)



- 1) GtxHu-001-D, Lot # 33-66-112613
- 2) Mol Weight Standard

Date: June 22, 2021 Signature: ______Quality Control

Phone: 919-831-2240 · www.ImmunoReagents.com · Fax: 919-831-2240