

6003 Chapel Hill Road · Suite # 153 · Raleigh, NC 27607

Certificate of Analysis

Product:	Donkey anti-Rabbit IgG (H&L), Affinity Pure
Part Number: Lot Number: Concentration: Amount: Form: Purification: Purify: Immunogen: Buffer: Preservative: Storage: Expiration Date:	DkxRb-003-D 80-86-042621 2.5 mg/ml (E ^{1%} at 280 nm = 13.0) 2.0 mg Clear, colorless liquid, 0.2 μ m filtered Affinity purified using solid phase Rabbit IgG \geq 95% based on SDS-PAGE Purified Rabbit IgG, whole molecule 10 mM Sodium Phosphate, 0.15 M Sodium Chloride, pH 7.2 0.05 % (w/v) Sodium Azide 2-8 ° C April 26, 2024
Specificity:	 Based on IEP, this antibody reacts with: heavy (γ) chains on Rabbit IgG light chains on all Rabbit immunoglobulins
Cross Reactivity:	Based on IEP, no reactivity is observed to:non-immunoglobulin Rabbit serum proteins
Country of Origin:	Donkey serum was obtained from healthy animals of US origin and under the care of a registered veterinarian.
Applications:	The optimal working dilution should be determined by the investigator. Suitable for use in: ELISA, IP, IDA, WB, IHC, Absorption
Disclaimer:	<i>For in vitro Laboratory Use Only.</i> Not for diagnostic or therapeutic use. Not for Human or animal consumption. The proper selection and use of our products is the sole responsibility of the end user, and therefore we can offer no guarantee to a specific experimental outcome. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license under any patent of ImmunoReagents, Inc.

Note: This product was produced by the dilution of part # DkxRb-003-D, lot # 67-97-090419

Product may not be resold or modified for resale without prior written approval of ImmunoReagents, Inc.

*ImmunoReagents,Inc.

6003 Chapel Hill Road · Suite # 153 · Raleigh, NC 27607

Testing Results

PRODUCT: PART NUMBER: LOT NUMBER:	Donkey anti-Rabbit IgG (H&L), Affinity Pure DkxRb-003-D 80-86-042621	
<u>TESTS</u>	SPECIFICATIONS	RESULTS
Appearance:	Clear and Colorless	PASS
Purity:	Two predominant bands at 50 and 24 kDa based on SDS PAGE	PASS
Specificity: Cross Reactivity:	 Based on IEP, this antibody has precipitin bands with: heavy (γ) chains on Rabbit IgG light chains on all Rabbit immunoglobulins Based on IEP, this antibody does not have precipitin bands to: non-immunoglobulin Rabbit serum proteins 	PASS PASS PASS
<u>SDS PAGE</u> 1 2	Immunoelectrophoresis IEP for Donkey anti-Rabbit IgG (H&L), Affinity Pure Part # DkxRb-003-D Lot # 67-97-09041	19
55kl	Da Normal Rabbit Ser	um
31kI	Da O Purified Rabbit I	gG
21kI	Da	
	03-D, Lot # 67-97-090419 ht Standard	
Date: April 27, 2021	Signature:Quality Control	