

6003 Chapel Hill Road · Suite # 153 · Raleigh, NC 27607

Certificate of Analysis

Product:	Donkey anti-Sheep IgG (H&L), Affinity Pure	
Part Number: Lot Number: Concentration: Amount: Form: Purification: Purity: Immunogen: Buffer: Preservative: Storage: Expiration Date:	DkxSh-003-D 69-159-042720 10.0 mg/ml (E ^{1%} at 280 nm = 13.0) 2.0 mg Clear, colorless liquid, 0.2 μ m filtered Affinity purified using solid phase Sheep IgG \geq 95% based on SDS-PAGE Purified Sheep IgG, whole molecule 10 mM Sodium Phosphate, 0.15 M Sodium Chloride, pH 7.2 0.05 % (w/v) Sodium Azide 2-8 ° C 04-27-2023	
Specificity:	 Based on IEP, this antibody reacts with: heavy (γ) chains on Sheep IgG light chains on all Sheep immunoglobulins 	
Cross Reactivity:	Based on IEP, no reactivity is observed to:non-immunoglobulin Sheep serum proteins	
Country of Origin:	Donkey serum was obtained from healthy animals of US origin and under the care of a registered veterinarian.	
Applications:	The optimal working dilution should be determined by the investigator. Suitable for use in: ELISA, IP, IDA, WB, IHC, Absorption	
Disclaimer:	<i>For in vitro Laboratory Use Only.</i> Not for diagnostic or therapeutic use. Not for Human or animal consumption. The proper selection and use of our products is the sole responsibility of the end user, and therefore we can offer no guarantee to a specific experimental outcome. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license under any patent of ImmunoReagents, Inc.	
	Product may not be resold or modified for resale without prior written approval of ImmunoReagents, Inc.	



6003 Chapel Hill Road · Suite # 153 · Raleigh, NC 27607

Testing Results

PRODUCT: PART NUMBER: LOT NUMBER:	Donkey anti- Sheep IgG (H&L), Affinity Pure DkxSh-003-D 69-159-042720	
<u>TESTS</u>	SPECIFICATIONS	<u>RESULTS</u>
Appearance:	Clear and Colorless	PASS
Purity:	Two predominant bands at 50 and 24 kDa based on SDS PAGE	PASS
Specificity:	 Based on IEP, this antibody has precipitin bands with: heavy (γ) chains on Sheep IgG light chains on all Sheep immunoglobulins 	PASS PASS
Cross Reactivity:	 Based on IEP, this antibody does not have precipitin bands to: non-immunoglobulin Sheep serum proteins 	PASS
SDS PAGE	<u>Immunoelectrophoresis</u>	
55k	Da IEP for Donkey anti-Sheep IgG (H&L), Affinity Da Part #: DkxSh-003-D Lot #: 69-159-	
	kDa Normal Sheep Ser	um
	DkxSh-C O Purified Sheep	Ĩ

- 1) DkxSh-003-D, Lot # 69-159-042720
- 2) Mol Weight Standard

Date: May 5, 2020

Signature: Quality Control